

Annexin A4

Catalog No: #AP80952



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Annexin A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	35-beta calcimedlin,Annexin IV,Annexin-4,Carbohydrate-binding protein p33/p41,Chromobindin-4,Endonexin I,Lipocortin IV,P32.5,PP4-X,Placental anticoagulant protein II,PAP-II,Protein II,ANX4
Accession No.	Uniprot ID:P09525
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.