

Transforming growth factor-beta-induced protein ig-h3

Catalog No: #AP80964



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor-beta-induced protein ig-h3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Beta ig-h3,MP70,MP78,BIGH3
Accession No.	Uniprot ID:P55906
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.