

BH3-interacting domain death agonist

Catalog No: #AP80970



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	BH3-interacting domain death agonist
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	p22 BID,BID
Accession No.	Uniprot ID:P55957
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.