

Sodium/potassium-transporting ATPase subunit alpha-3

Catalog No: #AP80982



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium/potassium-transporting ATPase subunit alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Na(+)/K(+) ATPase alpha-3 subunit,Na(+)/K(+) ATPase alpha(III) subunit,Sodium pump subunit alpha-3
Accession No.	Uniprot ID:Q6PIC6
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.