

B-cell receptor-associated protein 31

Catalog No: #AP80986



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	B-cell receptor-associated protein 31
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BCR-associated protein 31,p28,Bap31
Accession No.	Uniprot ID:Q61335
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.