

Cyclic AMP-dependent transcription factor ATF-4

Catalog No: #AP81002



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cyclic AMP-dependent transcription factor ATF-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	cAMP-dependent transcription factor ATF-4,Activating transcription factor 4,C/EBP-related ATF,C/ATF,Tax-responsive enhancer element-binding protein 67 homolog,TaxREB67 homolog
Accession No.	Uniprot ID:Q06507
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.