

Brain-specific serine protease 4

Catalog No: #AP81010



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Brain-specific serine protease 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	BSSP-4, Serine protease 22, Serine protease 26, Tryptase epsilon, Bssp4, Prss26
Accession No.	Uniprot ID:Q9ER10
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.