

Large proline-rich protein BAG6

Catalog No: #AP81044



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Large proline-rich protein BAG6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	BCL2-associated athanogene 6,HLA-B-associated transcript 3,Bat3
Accession No.	Uniprot ID:Q6MG49
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.