

Beta-secretase 1

Catalog No: #AP81068



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Beta-secretase 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Aspartyl protease 2,ASP2,Beta-site amyloid precursor protein cleaving enzyme 1,Beta-site APP cleaving enzyme 1,Memapsin-2,Membrane-associated aspartic protease 2,Bace |
| Accession No. | Uniprot ID:P56819 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.