

Voltage-dependent P/Q-type calcium channel subunit alpha-1A

Catalog No: #AP81474



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Voltage-dependent P/Q-type calcium channel subunit alpha-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Brain calcium channel I, BI, Calcium channel, L type, alpha-1 polypeptide isoform 4, Voltage-gated calcium channel subunit alpha Cav2.1, Caca1a, Cach4, Cacn3, Cacn1a4, Ccha1a
Accession No.	Uniprot ID:P97445
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.