

GST-Tag Monoclonal Antibody(1B10)

Catalog No: #44009



Package Size: #44009-1 50ul #44009-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GST-Tag Monoclonal Antibody(1B10)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1B10
Isotype	IgG1
Applications	WB
Specificity	The antibody detects C-terminal, internal, and N-terminal GST fusion proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	GST-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



0.5ug GST fusion protein+ Primary antibody dilution at

1) 1:5000

2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.