

# HSV-Tag Monoclonal Antibody(9D7)

Catalog No: #44016



Package Size: #44016-1 50ul #44016-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HSV-Tag Monoclonal Antibody(9D7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	9D7
Isotype	IgG1
Applications	WB
Specificity	The antibody detects HSV tag fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	HSV-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

## Images



1ug HSV fusion protein+ Primary antibody dilution at  
 1) 1:5000  
 2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.