

HA-Tag Monoclonal Antibody(1B10)

Catalog No: #44028



Package Size: #44028-1 50ul #44028-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | HA-Tag Monoclonal Antibody(1B10) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 1B10 |
| Isotype | IgG1 |
| Applications | WB |
| Specificity | The antibody detects C-terminal, internal, and N-terminal HA-tag fusion proteins. |
| Immunogen Type | Peptide |
| Immunogen Description | Synthetic Peptide |
| Target Name | HA-Tag |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of HA Recombinant protein, diluted at 1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.