

# Myc-Tag Monoclonal Antibody(3E8)

Catalog No: #44029



Package Size: #44029-1 50ul #44029-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Myc-Tag Monoclonal Antibody(3E8)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3.00E+08
Isotype	IgG1
Applications	WB
Specificity	The antibody can detect C-Myc fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Myc-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of c-Myc Recombinant protein, diluted at 1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.