

# CD1 Monoclonal Antibody(9H6)

Catalog No: #44036



Package Size: #44036-1 50ul #44036-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD1 Monoclonal Antibody(9H6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	9H6
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous CD1 proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	CD1
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Immunohistochemistry: 1:50-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.