

## GFP Monoclonal Antibody(Plant Specific)

Catalog No: #44050



Package Size: #44050-1 50ul #44050-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	GFP Monoclonal Antibody(Plant Specific)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Plant Specific
Isotype	IgG1
Applications	WB IP
Specificity	The antibody detects GFP, EGFP and GFP, EGFP tag fusion proteins in Plants.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	GFP
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

Immunoprecipitation: 1:200

## Images



1) Input: Hela Cell Lysate  
 2) IP product: IP dilute 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.