

# AmCyan Monoclonal Antibody(Mix)

Catalog No: #44065



Package Size: #44065-1 50ul #44065-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AmCyan Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB
Specificity	The antibody detects AmCyan & AmCyan Tag protein.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	AmCyan
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of recombinant AmCyan protein, diluted at 1:5000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.