

HSP90 α Monoclonal antibody(7D6)

Catalog No: #44078



Package Size: #44078-1 50ul #44078-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|-----------------------|---|
| Product Name | HSP90 α Monoclonal antibody(7D6) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 7D6 |
| Isotype | IgG1 |
| Applications | WB IHC |
| Species Reactivity | Hu Rt Ms |
| Specificity | The antibody detects endogenous HSP90 α protein |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant Protein |
| Target Name | HSP90 α |
| Other Names | HSP90 α |
| Calculated MW | 95kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details


Western blotting: 1:500-1:2000

Immunohistochemistry: 1:100-1:300

Images



Western blot analysis of 1)Hela, 2)Mouse Brain Tissue, 3) Rat Brain Tissue with HSP90 α Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma using HSP90 α Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.