

# PI3 Kinase P85α Monoclonal Antibody(3B7)

Catalog No: #44085



Package Size: #44085-1 50ul #44085-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PI3 Kinase P85α Monoclonal Antibody(3B7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3B7
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Rt Ms
Specificity	PI3 Kinase P85α Mouse Monoclonal Antibody(3B7) detects endogenous levels of PI3 Kinase P85α
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	PI3 Kinase P85α
Calculated MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:100-1:200


## Images



Western blot analysis of 1)3T3, 2) Rat LiverTissue with PI3 Kinase P85α Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat Liver Tissue using PI3 Kinase P85  $\alpha$  Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using PI3 Kinase P85  $\alpha$  Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.