

## TBP/TATA Binding Protein Monoclonal Antibody(1F6)

Catalog No: #44097



Package Size: #44097-1 50ul #44097-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TBP/TATA Binding Protein Monoclonal Antibody(1F6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1F6
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	TBP/TATA Binding Protein Mouse Monoclonal Antibody(1F6) detects endogenous levels of TBP/TATA Binding Protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	TBP/TATA Binding Protein
Calculated MW	35-45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using TBP/TATA Binding Protein mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using TBP/TATA Binding Protein Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.