

Vesicular integral-membrane protein VIP36

Catalog No: #AP84907



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vesicular integral-membrane protein VIP36
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lectin mannose-binding 2, Vesicular integral-membrane protein 36, VIP36
Accession No.	Uniprot ID:Q9DBH5
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.