

Transferrin Monoclonal Antibody(10G3)

Catalog No: #44109



Package Size: #44109-1 50ul #44109-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|-----------------------|--|
| Product Name | Transferrin Monoclonal Antibody(10G3) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 10G3 |
| Isotype | IgG1 |
| Applications | WB IHC |
| Species Reactivity | Hu |
| Specificity | Transferrin Mouse Monoclonal AntibodyB£B"10G3B£B© detects endogenous levels of Transferrin |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant Protein |
| Target Name | Transferrin |
| Calculated MW | 77kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western blotting: 1:2000-1:5000

Immunohistochemistry: 1:100-1:200

Images



Western blot analysis of Human Serum with Transferrin Mouse mAb diluted at 1:2,000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.