

Transferrin Monoclonal Antibody(3A5)

Catalog No: #44110



Package Size: #44110-1 50ul #44110-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transferrin Monoclonal Antibody(3A5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3A5
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB£B™3A5B£B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Transferrin
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:100-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.