

# Transferrin Monoclonal Antibody(1B12)

Catalog No: #44112



Package Size: #44112-1 50ul #44112-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transferrin Monoclonal Antibody(1B12)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1B12
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB&B™1B12B&B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Transferrin
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Immunohistochemistry: 1:100-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.