

# HIF-1 $\beta$ /ARNT Monoclonal Antibody(4C5)

Catalog No: #44119



Package Size: #44119-1 50ul #44119-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	HIF-1 $\beta$ /ARNT Monoclonal Antibody(4C5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4C5
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Ms
Specificity	HIF-1 $\beta$ /ARNT Mouse Monoclonal Antibody(4C5) detects endogenous levels of HIF-1 $\beta$ /ARNT
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	HIF-1 $\beta$ /ARNT
Calculated MW	87kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:100-1:200

## Images



Western blot analysis of Mouse Brain Tissue with HIF-1  $\beta$ /ARNT Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Mouse BrainTissue using HIF-1  $\beta$ /ARNT Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.