

IL-8 Mouse Monoclonal Antibody(13F8)

Catalog No: #44138



Package Size: #44138-1 50ul #44138-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL-8 Mouse Monoclonal Antibody(13F8)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	13F8
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	IL-8 Mouse Monoclonal Antibody(13F8) detects endogenous levels of IL8
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	IL-8
Other Names	IL8
Calculated MW	11kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using IL-8 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using IL-8 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.