

## Cyclin B1 Mouse Monoclonal Antibody(1A5)

Catalog No: #44177



Package Size: #44177-1 50ul #44177-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Cyclin B1 Mouse Monoclonal Antibody(1A5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1A5
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	Cyclin B1 Mouse Monoclonal Antibody(1A5) detects endogenous levels of CCNB1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	Cyclin B1
Other Names	CCNB1
Calculated MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-1:2000

Immunohistochemistry: 1:100-1:200

## Images



Western blot analysis of 1) Hela Cell Lysate, 2) 3T3 Cell Lysate, 3) Rat Brain Tissue Lysate using LC3A Mouse mAb diluted at 1:1000.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Cyclin B1 Mouse mAb diluted at 1:200.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.