

Neutrophil gelatinase-associated lipocalin

Catalog No: #AP85722



Package Size: #AP85722-1 50ug #AP85722-2 100ug #AP85722-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neutrophil gelatinase-associated lipocalin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NGAL,Lipocalin-2,SV-40-induced 24P3 protein,Siderocalin LCN2,p25
Accession No.	Uniprot ID:P11672
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.