

Collagen IV Mouse Monoclonal Antibody(8E5)

Catalog No: #44258



Package Size: #44258-1 50ul #44258-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Collagen IV Mouse Monoclonal Antibody(8E5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8.00E+05
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Collagen IV protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Collagen IV
Target Name	Collagen IV
Other Names	Collagen alpha-1(IV) chain [Cleaved into: Arresten]
Calculated MW	161kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using Collagen IV Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using Collagen IV Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.