

## PAK4/5/6 (Phospho-Ser474) Antibody

Catalog No: #12417

Package Size: #12417-1 50ul #12417-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PAK4/5/6 (Phospho-Ser474) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PAK4/5/6 (Phospho-Ser474) Antibody detects endogenous levels of PAK4/5/6 only when phosphorylated at Ser474
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PAK4/5/6 (Phospho-Ser474)
Target Name	PAK4/5/6
Modification	Phospho
Other Names	PAK4, KIAA1142, p21(CDKN1A)-activated kinase 4, PAK-4, p21-activated kinase 4
Accession No.	Swiss-Prot#: O96013NCBI Gene ID: 10298
Target Species	human
Calculated MW	64kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

## Images



Western Blot analysis of cell lysates from HepG2 and HepG2 UV treated, using PAK4/5/6 (Phospho-Ser474) Antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.