

p27 Kip1 (Phospho-Ser178) Antibody

Catalog No: #12430



Package Size: #12430-1 50ul #12430-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	p27 Kip1 (Phospho-Ser178) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	p27 Kip1 (Phospho-Ser178) Antibody detects endogenous levels of p27 Kip1 only when phosphorylated at Ser178
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human p27 Kip1 (Phospho-Ser178)
Target Name	p27 Kip1
Modification	Phospho
Other Names	CDKN1B, KIP1, MEN1B, P27 Kip1, p27KIP1, MEN4, CDKN4
Accession No.	Swiss-Prot#: P46527NCBI Gene ID: 1027
Target Species	human
Calculated MW	27kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.