

MAX (Phospho-Ser2) Antibody

Catalog No: #12449



Package Size: #12449-1 50ul #12449-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAX (Phospho-Ser2) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	MAX (Phospho-Ser2) Antibody detects endogenous levels of MAX only when phosphorylated at Ser2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MAX (Phospho-Ser2)
Target Name	MAX
Modification	Phospho
Other Names	MAX, BHLHd4, BHLHd5, BHLHd6, BHLHd8, MYC associated factor X, ORF1, Protein max, BHLHd7, MAX protein, Myc-associated factor X
Accession No.	Swiss-Prot#: P61244NCBI Gene ID: 4149
Target Species	human
Calculated MW	34kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.