

Rel (Phospho-Ser492/460) Antibody

Catalog No: #12484



Package Size: #12484-1 50ul #12484-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rel (Phospho-Ser492/460) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	Rel (Phospho-Ser492/460) Antibody detects endogenous levels of Rel only when phosphorylated at Ser492/460
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Rel (Phospho-Ser492/460)
Target Name	Rel
Modification	Phospho
Other Names	REL, C-Rel, C-Rel proto-oncogene protein, I-Rel, Proto-oncogene c-Rel
Accession No.	Swiss-Prot#: Q04864NCBI Gene ID: 5966
Target Species	human
Calculated MW	68kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.