

ABL1/2 (Phospho-Tyr89/Tyr70) Antibody

Catalog No: #12488



Package Size: #12488-1 50ul #12488-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABL1/2 (Phospho-Tyr89/Tyr70) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ABL1/2 (Phospho-Tyr89/Tyr70) Antibody detects endogenous levels of ABL1/2 only when phosphorylated at Tyr89/Tyr70
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ABL1/2 (Phospho-Tyr89/Tyr70)
Target Name	ABL1/2
Modification	Phospho
Other Names	ABL1, Bcr/c-abl oncogene protein, Bcr-Abl, Bcr/abl, C-ABL, ABL, JTK7, p150, Tyrosine-protein kinase ABL1, Proto-oncogene c-Abl, V-abl
Accession No.	Swiss-Prot#: P00519/P42684NCBI Gene ID: 25/27
Target Species	human
Calculated MW	135kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.