

FGFR2 (Phospho-Tyr769) Antibody

Catalog No: #12503

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FGFR2 (Phospho-Tyr769) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen
Applications	WB
Species Reactivity	Hu Ms
Specificity	FGFR2 (Phospho-Tyr769) Antibody detects endogenous levels of FGFR2 only when phosphorylated at Tyr769
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FGFR2 (Phospho-Tyr769)
Target Name	FGFR2
Modification	Phospho
Other Names	FGFR2, BEK, BBDS, Bacteria-expressed kinase, CD332 antigen, CFD1, Craniofacial dysostosis 1, CEK3, Fgf bek, FGF receptor, Hydroxyaryl-protein kinase, Jackson-Weiss syndrome, JWS, FGF Receptor 2, KSAM, KGFR, Soluble FGFR4 variant 4, BFR-1, TK14, TK25, ...
Accession No.	Swiss-Prot#: P21802NCBI Gene ID: 2263
Target Species	human
Calculated MW	92kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.