

FRS2 (Phospho-Tyr196) Antibody

Catalog No: #12539

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FRS2 (Phospho-Tyr196) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	FRS2 (Phospho-Tyr196) Antibody detects endogenous levels of FRS2 only when phosphorylated at Tyr196
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FRS2 (Phospho-Tyr196)
Target Name	FRS2
Modification	Phospho
Other Names	FRS2, FGFR substrate 2, FRS2A, FGFR signalling adaptor, FGFR-signaling adaptor SNT, FRS2alpha, SNT-1, SNT1, SNT
Accession No.	Swiss-Prot#: Q8WU20NCBI Gene ID: 10818
Target Species	human
Calculated MW	57kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.