

C-RAF (Phospho-Tyr340) Antibody

Catalog No: #12559

Package Size: #12559-1 50ul #12559-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-RAF (Phospho-Tyr340) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	C-RAF (Phospho-Tyr340) Antibody detects endogenous levels of C-RAF only when phosphorylated at Tyr340
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human C-RAF (Phospho-Tyr340)
Target Name	C-RAF
Modification	Phospho
Other Names	RAF1, C-Raf, C-raf-1, CRAF, NS5, Raf-1, Proto-oncogene c-RAF, Oncogene RAF1, RAF, Raf1 kinase
Accession No.	Swiss-Prot#: P04049NCBI Gene ID: 5894
Target Species	human
Calculated MW	73kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.