

# PDGFRa (Phospho-Tyr1018) Antibody

Catalog No: #12560

Package Size: #12560-1 50ul #12560-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PDGFRa (Phospho-Tyr1018) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PDGFRa (Phospho-Tyr1018) Antibody detects endogenous levels of PDGFRa only when phosphorylated at Tyr1018
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PDGFRa (Phospho-Tyr1018)
Target Name	PDGFRa
Modification	Phospho
Other Names	PDGFRA, CD140A, CD140a antigen, PDGFalphaR, PDGFR2, PDGFRalpha, AlphaPDGFR, RHEPDGFRA, PDGF receptor alpha, PDGFR-2, PDGFRA/BCR fusion, APDGFR, Pdgf alpha-receptor, PDGF-R-alpha, PDGFR-alpha
Accession No.	Swiss-Prot#: P16234NCBI Gene ID: 5156
Target Species	human
Calculated MW	140kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.