

BIK (Phospho-Ser35) Antibody

Catalog No: #12578



Package Size: #12578-1 50ul #12578-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BIK (Phospho-Ser35) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	BIK (Phospho-Ser35) Antibody detects endogenous levels of BIK only when phosphorylated at Ser35
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human BIK (Phospho-Ser35)
Target Name	BIK
Modification	Phospho
Other Names	BIK, Apoptosis inducer NBK, Bcl-2-interacting killer, BP4, Apoptosis-inducing NBK, BIP1, NBK
Accession No.	Swiss-Prot#: Q13323NCBI Gene ID: 638
Target Species	human
Calculated MW	30kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.