

## Fos (Phospho-Ser374) Antibody

Catalog No: #12599



Package Size: #12599-1 50ul #12599-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fos (Phospho-Ser374) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Fos (Phospho-Ser374) Antibody detects endogenous levels of Fos only when phosphorylated at Ser374
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Fos (Phospho-Ser374)
Target Name	Fos
Modification	Phospho
Other Names	FOS, Activator protein 1, AP-1, Cellular oncogene c-fos, Cellular oncogene fos, Proto-oncogene c-Fos, v-fos, C-FOS, G0s7, p55
Accession No.	Swiss-Prot#: P01100NCBI Gene ID: 2353
Target Species	human
Calculated MW	50kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.