

V-type proton ATPase catalytic subunit A

Catalog No: #AP88807



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	V-type proton ATPase catalytic subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	V-ATPase subunit A,V-ATPase 69 kDa subunit,Vacuolar proton pump subunit alpha,ATP6A1,ATP6V1A1
Accession No.	Uniprot ID:Q29048
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.