

NACA (Phospho-Thr159) Antibody

Catalog No: #12618



Package Size: #12618-1 50ul #12618-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NACA (Phospho-Thr159) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	NACA (Phospho-Thr159) Antibody detects endogenous levels of NACA only when phosphorylated at Thr159
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NACA (Phospho-Thr159)
Target Name	NACA
Modification	Phospho
Other Names	NACA, Allergen Hom s 2, Alpha-NAC, HSD48, NAC-alpha, NACA1
Accession No.	Swiss-Prot#: Q13765/Q9H009NCBI Gene ID: 4666/342538
Target Species	human
Calculated MW	24kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.