

# Vasopressin V1b receptor

Catalog No: #AP88908



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vasopressin V1b receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	V1bR,AVPR V1b,AVPR V3,Antidiuretic hormone receptor 1b,Vasopressin V3 receptor
Accession No.	Uniprot ID:P48974
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.