

PTP1B (Phospho-Tyr66) Antibody

Catalog No: #12630



Package Size: #12630-1 50ul #12630-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PTP1B (Phospho-Tyr66) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PTP1B (Phospho-Tyr66) Antibody detects endogenous levels of PTP1B only when phosphorylated at Tyr66
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PTP1B (Phospho-Tyr66)
Target Name	PTP1B
Modification	Phospho
Other Names	PTPN1, PTP1B, PTP-1B, PTPase 1B
Accession No.	Swiss-Prot#: P18031NCBI Gene ID: 5770
Target Species	human
Calculated MW	49kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western blot analysis PTP1B (Phospho-Tyr66) using MDA-MB-435 whole cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.