

Y-box-binding protein 3

Catalog No: #AP89002



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Y-box-binding protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cold shock domain-containing protein A,DNA-binding protein A,Single-strand DNA-binding protein NF-GMB,CSDA,DBPA
Accession No.	Uniprot ID:P16989
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.