

Vasoactive intestinal polypeptide receptor 2

Catalog No: #AP89004



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vasoactive intestinal polypeptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	VIP-R-2,Pituitary adenylate cyclase-activating polypeptide type III receptor,PACAP type III receptor
Accession No.	Uniprot ID:P35000
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.