

SRF (Phospho-Thr159) Antibody

Catalog No: #12642



Package Size: #12642-1 50ul #12642-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SRF (Phospho-Thr159) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	SRF (Phospho-Thr159) Antibody detects endogenous levels of SRF only when phosphorylated at Thr159
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SRF (Phospho-Thr159)
Target Name	SRF
Modification	Phospho
Other Names	SRF, Serum response factor, MCM1
Accession No.	Swiss-Prot#: P11831NCBI Gene ID: 6722
Target Species	human
Calculated MW	55kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.