

Vitamin D3 receptor

Catalog No: #AP89102



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vitamin D3 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	VDR,1,25-dihydroxyvitamin D3 receptor,Nuclear receptor subfamily 1 group I member 1,NR111
Accession No.	Uniprot ID:A3RGC1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.