

VHL (Phospho-Ser68) Antibody

Catalog No: #12654



Package Size: #12654-1 50ul #12654-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VHL (Phospho-Ser68) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	VHL (Phospho-Ser68) Antibody detects endogenous levels of VHL only when phosphorylated at Ser68
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human VHL (Phospho-Ser68)
Target Name	VHL
Modification	Phospho
Other Names	VHL, Elongin binding protein, HRCA1, RCA1, VHL1, Von Hippel-Lindau syndrome, Protein G7, PVHL
Accession No.	Swiss-Prot#: P40337NCBI Gene ID: 7428
Target Species	human
Calculated MW	24kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.