

53BP1 (Phospho-Ser1778) Antibody

Catalog No: #12657



Package Size: #12657-1 50ul #12657-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	53BP1 (Phospho-Ser1778) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	53BP1 (Phospho-Ser1778) Antibody detects endogenous levels of 53BP1 only when phosphorylated at Ser1778
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human 53BP1 (Phospho-Ser1778)
Target Name	53BP1
Modification	Phospho
Other Names	TP53BP1, 53BP1, p202, p53-binding protein 1, p53-binding protein 53bp1, p53BP1
Accession No.	Swiss-Prot#: Q12888NCBI Gene ID: 7158
Target Species	human
Calculated MW	214kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.